

BRCAaccuTest & BRCAaccuTest Plus

Genetic testing for BRCA1/2 mutations

For breast and ovarian cancer patients
For HBOC (Hereditary breast and ovarian cancer syndrome) patients
For family history or the age of breast cancer incidence


BRCAaccuTest™ is a reagent for producing libraries for analyzing the BRCA1/2 genes using the NGS (next-generation sequencing) method, which analyzes genomic DNA derived from blood or FFPE tissues.



Product Specification

	BRCAaccuTest™	BRCAaccuTest™ PLUS
certification	CE-IVD, KFDA Class III	CE-IVD
Technique	NGS (Next Generation Sequencing)	
Compatible platforms	Illumina / MiSeq, MiSeqDx	
Target Enrichment	Amplicon method	
Specimen	Blood	Blood and FFPE tissue, tissue biopsy
Quantity	24 test (1 test for control DNA)	6+6 test (1 test each for control DNA)
Target region	19.8 kb (all protein coding region + intron boundary 10 bp)	
Amplicons(primer pairs)	131 (2 pools)	146 (2 pools)
Mean amplicon size	234 bp	207 bp
Hands on time	4.5 hrs	5 hrs

Working Process




Product Performance


	Germline samples	Somatic samples
Compatible kit	BRCAaccuTest™ / BRCAaccuTest™ PLUS	BRCAaccuTest™ PLUS
Target coverage	Mean coverage 200X Minimum coverage 20X (MiSeq Reagent Nano kit v2, 300cycles / 24 tests)	Mean coverage 1000X Minimum coverage 200X (Using MiSeq Reagent Nano kit v2, 300cycles / 6 tests)
Uniformity	100% detected >0.2x mean coverage	100% detected >0.2x mean coverage
Accuracy	> 99.99%	> 99.99%
Target specificity	> 95% on target read counts	> 95% on target read counts
Sensitivity	Input gDNA 0.5 ng	VAF 5%

BRCAaccuTest™ & BRCAaccuTest™ PLUS Comparison

BRCAaccuTest™




BRCAaccuTest™ PLUS



BRCAaccuTest™ Validation Data

BRCAaccuTest™ gave the positive, negative, and overall diagnostic consistency over 99 % with Sanger sequencing at a significance level of 5%

Clinical Study Design



Diagnostic agreement : 100%

		Sanger Sequencing	
		P	N
BRCAaccuTest	P	TP	FP
	N	FN	TN

* True Positive (TP), False Positive (FP)
* False Negative (FN), True Negative (TN)

BRCA1 (BRCA2)		Sanger Sequencing		
		P	N	Sum
BRCAaccuTest	P	103	0	103
	N	0	103	103
Sum		103	103	206

* Test results in BRCA2 were identical to those in BRCA1
* Tested number of samples

GENE	PPA (95% CI)	NPA (95% CI)	OA (95% CI)
BRCA1	100 (0.9648, 1)	100 (0.9648, 1)	100 (0.9823, 1)
BRCA2	100 (0.9648, 1)	100 (0.9648, 1)	100 (0.9823, 1)

* Positive Percent Agreement (PPA) = $[\text{TP}/(\text{TP}+\text{FN})] \times 100$
* Negative Percent Agreement (NPA) = $[\text{TN}/(\text{TN}+\text{FP})] \times 100$
* Overall Agreement (OA) = $[(\text{TP}+\text{TN})/(\text{TP}+\text{FP}+\text{FN}+\text{TN})] \times 100$

Variants agreement : 100%

BRCA1 (BRCA2)		Sanger Sequencing			
		SNV	Duplication	Deletion	INDEL
BRCAaccuTest	SNV	1,629	0	0	0
	Duplication	0	12	0	0
	Deletion	0	0	49	0
	INDEL	0	0	0	14

* Number of detected variants

NGeneAnalySys™ Performance

Test Data		Test Number (Variants tested)		Concordance
		Manual Analysis	NGeneAnalySys™	
Analytical Validation	Cell lines	277 (219)	277 (219)	100%
	Clinical Samples	80 (78)	80 (78)	100%
Clinical Validation	Clinical Samples	206 (1,704)	206 (1,704)	100%
Sum		563 (2,001)	563 (2,001)	



Ordering Information

Cat no.	Products	Quantity (T)	Storage
NGB111VKA	<p>BRCAaccuTest™ (Approved CE-IVD and KFDA Class III)</p> <ul style="list-style-type: none"> - BRCA Amplicon Primer Mix 1 - BRCA Amplicon Primer Mix 2 - Hot Start PCR I Master Mix - Ligation Mix - Amplicon DNA Adapter - PCR II Barcoded Primer Mix N1~N24 - PCR II Mater Mix - Nuclease-Free Water - Control DNA - Resuspension Buffer - Cleanup Beads 	24	<p>-20°C (box package)</p> <p>4°C (separated component)</p>
NGB112VIA	<p>BRCAaccuTest™ PLUS (Approved CE-IVD)</p> <ul style="list-style-type: none"> - BRCA Primer Mix 1 - BRCA Primer Mix 2 - Hot Start PCR Polymerase - 2X PCR Buffer - End Repair Enzyme mix - End Repair Reaction Buffer - Ligation Master Mix - Ligation Enhancer - Adapter (for Illumina) - Uracil-DNA glycosylase (UDG) - PCR II Mater Mix - Barcoded Primer (Index1~Index12) - Universal Primer - Control DNA 	12	<p>-20°C (box package)</p>
NGB311VWW	<p>NGeneAnalySys™ v.1.1 (Software for only germline analysis)</p>		
NGB312VWW	<p>NGeneAnalySys™ v.1.3 (Software for both germline and somatic analysis)</p>		

This product can be used for *in vitro* diagnostic purpose in authorized countries

All products of NGeneBio Co., Ltd. are manufactured under a strict quality control process. NGeneBio Co., Ltd. warrants the quality of the product during the warranty period (marked on the product).